

# SHUR/Mark<sup>®</sup> Cassette Sleeves

*Prepackaged Cassette Sleeves insert directly  
into SHUR/Mark<sup>®</sup> Hoppers*




**TBS<sup>®</sup>**

# How To Load SHUR/Mark® Cassette Sleeves



Specialty dispensing hopper (SM-TSDH) accepts prepkg'd cassette sleeves.




Insert cassette sleeves into hopper sitting on 4-position hopper rack (SM-HR).



Place loaded hopper on SHUR/Mark® and pull tape down to seat cassettes.

## Ordering Information

	White	Blue	Green	Gray	Pink	Peach	Yellow	Tan	Aqua	Lilac
<b>Cat #</b>	<b>U/M</b>		<b>Description</b>							
SM-PCTS-(color)	cs		<b>SHUR/Mark®</b> Process & Embedding cassettes. Prepackaged <b>SHUR/Track™</b> taped cassette sleeves for direct insertion into SM-TSDH dispensing hopper receptacle. Colors Available: White ( <b>W</b> ), Blue ( <b>B</b> ), Green ( <b>G</b> ), Gray ( <b>GY</b> ), Pink ( <b>P</b> ), Peach ( <b>PE</b> ), Yellow ( <b>Y</b> ), Tan ( <b>T</b> ), Aqua ( <b>A</b> ) & Lilac ( <b>L</b> ); Packaged (20 sleeves of 50 cassette bottoms & 1 box of 1000 tops)/cs							
SM-BCTS-(color)	cs		<b>SHUR/Mark®</b> Biopsy cassettes. Prepackaged <b>SHUR/Track™</b> biopsy taped cassette sleeves for direct insertion into SM-TSDH dispensing hopper receptacle. Colors Available: White ( <b>W</b> ), Blue ( <b>B</b> ), Green ( <b>G</b> ), Gray ( <b>GY</b> ), Pink ( <b>P</b> ), Peach ( <b>PE</b> ), Yellow ( <b>Y</b> ), Tan ( <b>T</b> ), Aqua ( <b>A</b> ) & Lilac ( <b>L</b> ); Packaged (20 sleeves of 50 cassette bottoms & 1 box of 1000 tops)/cs							
SM-TSDH	ea		Dispensing hopper for <b>SHUR/Track™</b> Cassette Tape Sleeves, reusable metal, accepts prepackaged taped cassette sleeves SM-PCTS-(color) and SM-BCTS-(color) biopsy cassettes, 50 cassette capacity							
SM-HR	ea		Hopper Rack, 4 position for holding additional (SM-TSDH) hoppers to facilitate loading of prepackaged cassette sleeves (SM-PCTS-(color)), hardware included for table or wall mounting							



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

[www.trianglebiomedical.com](http://www.trianglebiomedical.com)

Distributed By:

**TBS® TRIANGLE BIOMEDICAL SCIENCES, INC.**

3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX:919.384.9595 • EMAIL: [tbs@trianglebiomedical.com](mailto:tbs@trianglebiomedical.com)

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability.

SHUR/Track is a trademark, and SHUR/Mark® & TBS® are registered trademarks of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.

Litho in U.S.A. Printed September 2006

Hopper Patent #6.098.839, Patent Applied for Tape Mechanism

Bulletin 23.2.3