



Date Issued: **February 3, 2011**

Supersedes Date: **November 20, 2008**

## MATERIAL SAFETY DATA SHEET

IDENTIFICATION AND EMERGENCY INFORMATION				
<b>Product Name</b>	DabOut			
<b>Product Code (Catalog Number)</b>	DO-C, DO-LVS-C			
<b>Product Category (Chemical Name)</b>	SolvCycle A			
<b>Product Appearance and Odor</b>	Black granular or powder, no odor.			
<b>Emergency Notification</b>	Chem-Tel 1-800-255-3924			
COMPONENTS AND HAZARD INFORMATION				
<b>Component</b>	<b>CAS No.</b>	<b>%</b>	<b>OSHA TWA</b>	<b>OSHA STEL</b>
Activated Carbon	7440-44-0		NA	NA
PHYSICAL DATA				
<b>THE FOLLOWING DATA ARE APPROXIMATE OR TYPICAL VALUES AND SHOULD NOT BE USED FOR PRECISE DESIGN PURPOSES.</b>				
<b>Boiling Range</b> NA	<b>Specific Gravity</b> 1.9-2.2	<b>Molecular Weight</b> NL	<b>Viscosity</b> NA	
<b>Percent Volatile by Volume</b> NA	<b>Vapor Pressure @20 C degrees</b> NA	<b>Vapor Density (Air = 1)</b> NA	<b>Pour, Congealing or Melting Point</b> NA	
<b>Solubility in Water @1 ATM and 25°C (77°F)</b> NA		<b>Evaporation Rate @1 ATM and 25°C (77°F) (n-Butyl Acetate =1)</b> NA		<b>pH</b> NA
FIRE AND EXPLOSION HAZARD INFORMATION				
<b>Flash Point (Minimum)</b> NA	<b>Autoignition Temperature</b> 350 °C	<b>Test Mode</b> NL	<b>LEL</b> NL	<b>UEL</b> NL
<b>Handling Precautions</b> None, as long as the cartridge is in tact. Dust from broken cartridges may cause the following: Eye irritation, minor irritation to respiratory system. (ACGIH, OSHA and other TLV are not applicable for this product)				
<b>Flammable or Explosive Limits (Approximate Percent by Volume in Air)</b> Stable under normal handling and storage conditions.				

NA = not applicable

NL = not listed

NE = not established

NNR = not normally required

<p><b>Extinguishing Media</b> Flood with plenty of water or inert gas such as N<sub>2</sub> and CO<sub>2</sub>.</p>
<p><b>Fire Fighting Procedures</b> None. Over heating can produce semi-liquid melt. Molten material can also contribute to spread of fire.</p>
<p><b>Decomposition Products Under Fire Condition</b> Possible fire combustion products: Alkybenzenes, vinylbenzenes, naphthalene, benzaldehydes, phenol, carbon dioxide, water, organic amines, ammonia, chlorine, nitrogen oxide, methyl chloride and minor amounts of hydrogen cyanide, or aldehydes.</p>
<p><b>“Empty” Container Warning</b> NE</p>
<p><b>Unusual Fire and Explosion Hazards</b> Warning: Strong oxidizing agents such as nitric acid or ozone may, under certain conditions, degrade the contents of the cartridge up to a explosive reaction.</p>
<p><b>HEALTH AND HAZARD INFORMATION</b></p>
<p><b>Variability Among Individuals</b> NE</p>
<p><b>Signs and symptoms of overexposure</b> Slight irritation of eyes and nose may result from contact.</p>
<p><b>Effects of Overexposure</b> NA</p>
<p><b>Nature of Hazard and Toxicity Information</b> NNR</p>
<p><b>Pre-existing Medical Conditions which may be Aggravated by Exposure</b> NA</p>
<p><b>Threshold Limit Value</b> NE</p>
<p><b>Health Hazards</b> Eye: May cause severe eye irritation. May cause moderate corneal injury. Effects are likely to heal Skin: Prolonged or repeated exposure not likely to cause significant skin irritation. May cause more severe response if skin is abraded. Ingestion: vapors are unlikely due to physical properties.</p>
<p><b>PRIMARY ROUTES OF ENTRY AND EMERGENCY AND FIRST AID PROCEDURES</b></p>
<p><b>Eye Contact</b> Flush with plenty of water for as least 15 minutes.</p>
<p><b>Skin</b> Wash off with soap and water.</p>
<p><b>Inhalation</b> No adverse effects anticipated by this route of exposure.</p>

NA = not applicable

NL = not listed

NE = not established

NNR = not normally required

**Ingestion**

No adverse effects anticipated by this route of exposure.

Note to physician: No specific antidote. Supportive care. Treatment based on the judgement of the physician in response to reaction of the patient.

**SEEK MEDICAL TREATMENT IF DISCOMFORT PERSISTS.**

**REACTIVITY**

<b>Stability</b> Stable under normal conditions	<b>Incompatibility</b> Strong oxidizing agents such as nitric acid, ketone or ozone.	<b>Hazardous Decomposition Products</b> CO
<b>Conditions to Avoid</b> Generation of dust during and fines during handling of a broken cartridge.	<b>Hazardous Polymerization</b> Will not occur	<b>Materials to Avoid</b> See incompatibility statement

**ENVIRONMENTAL INFORMATION****Steps to be Taken in Case Material is Released or Spilled**

Wear protective breathing filter and sweep up.

**Disposal**

Dispose in accordance with local, state and federal regulations. Incineration.

**Protective Gear**

Rubber gloves, Protective clothing, Safety glasses, If dust from a cartridge is present use a NIOSH approved particulate filter. Good general ventilation.

**Precautions to be taken in handling and storage**

Avoid breaking open the cartridge as dust may occur. Avoid possible sources of ignition.

**THE FOLLOWING INFORMATION MAY BE USEFUL IN COMPLYING WITH VARIOUS STATE AND FEDERAL LAWS AND REGULATIONS UNDER VARIOUS ENVIRONMENTAL STATUTES.**

**Reportable Quantity (RQ), EPA Regulation 40 CFR 302 (CERCLA Section 102)**

NL

**Threshold Planning Quantity (TPQ), EPA Regulation 40 CFR 355 (SARA Sections 301-304)**

NL

**Toxic Chemical Release Reporting, EPA Regulation 40 CFR 372 (SARA Section 313)**

NL

**Hazardous Chemical Reporting, EPA Regulation 40 CFR 370 (SARA Sections 311-312)**

NL

**PROTECTION AND PRECAUTIONS****Ventilation**

Good general ventilation.

**Respiratory Protection**

No respiratory protection should be required.

NA = not applicable

NL = not listed

NE = not established

NNR = not normally required

<b>Protective Gloves</b> Not needed with normal use.
<b>Eye Protection</b> Optional
<b>Other Protective Equipment</b> Safety shower and eye bath.
<b>TRANSPORTATION AND OSHA RELATED LABEL INFORMATION</b>
<b>Transportation Incident Information</b> Notify EPA.
<b>U.S. DOT Hazardous Materials Shipping Description</b> NL
<b>OSHA Required Label Information</b> NL

**TRIANGLE BIOMEDICAL SCIENCES (TBS) provides the information contained herein in good faith, but makes no representation as to its comprehensiveness or accuracy. Individuals receiving this information must exercise their independent judgement in determining its appropriateness for a particular purpose. TBS makes no representation or warranties, either expressed or implied of merchantability, fitness for particular purposes with respect to the information set forth herein or to which the information refers. Accordingly, TBS will not be responsible for damages resulting from the use of or reliance upon this information.**