

# Forceps Warmer



**TBS**<sup>®</sup>

Has Your

# Forceps Warmer


Seen Better Days...



*If so, try TBS' indestructible, all metal, space saving unit. Three heated wells prevent wax build-up and allow easy tissue orientation at embedding centers.*

## Ordering Information:

Cat #	U/M	Description
FW-120	ea	Forceps Warmer, w/ 3 heated wells; 100-120 VAC, 50/60Hz, 0.5A; 5"W x 3 1/4"D x 3"H
FW-220	ea	Forceps Warmer, w/ 3 heated wells; 220-240 VAC, 50/60Hz, 0.25A; 5"W x 3 1/4"D x 3"H



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

[www.trianglebiomedical.com](http://www.trianglebiomedical.com)

Distributed By:

**TBS** TRIANGLE BIOMEDICAL SCIENCES, INC.

3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX:919.384.9595 • EMAIL: [tbs@trianglebiomedical.com](mailto:tbs@trianglebiomedical.com)

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. TBS is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.

Litho in U.S.A. © February 2008

Bulletin # 13..3.2