

Paraffin Dispenser



6.25 Gallons of Liquid Paraffin on Demand



TBS Anticipates Your Needs . . .

6.25 Gallons of Liquid Paraffin on Demand.

Large Capacity

6.25 gallon reservoir to meet your liquid paraffin requirements

Convenience

High spigot positioning accomodates a variety of large vessels

Microprocessor Controlled

Accurate temperature control and uniformity for today's paraffins

Clear Lexan Lid

Permits easy monitoring and refilling of paraffin level


Security

Advanced software & mechanical devices assure unparalled safety



Ordering Information:

Cat #	U/M	Description
PD-120	ea	Paraffin Dispenser, microprocessor controlled; 6.25 gallon capacity; 100-120VAC, 50/60Hz, 5A; ambient to 65°C; 14"W x 23"H x 15"D
PD-220	ea	Paraffin Dispenser, microprocessor controlled; 6.25 gallon capacity; 220-240VAC, 50/60Hz, 2.5A; ambient to 65°C; 14"W x 23"H x 15"D



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

www.trianglebiomedical.com

Distributed By:

TBS TRIANGLE BIOMEDICAL SCIENCES, INC.

3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX:919.384.9595 • EMAIL: tbs@trianglebiomedical.com

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. TBS is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.

Litho in U.S.A. © February 2006

Bulletin # 10.1.7