

# Slide Dryer II



*Makes Slide Drying a Breeze!*

**TBS**<sup>®</sup>

# Slide Dryer II..


*The Convenient Way To Dry Slides.*



- Temperature control to 75° C
- Lexan front for easy viewing
- Thermal switch for overheat protection
- Accommodates two slide stainer racks

## Ordering Information:

Cat #	U/M	Description
SD-II-120	ea	Slide Dryer II, forced air; 100-120VAC, 50/60Hz, 2A; ambient to 75°C; 11"W x 10"D x 5"H
SD-II-220	ea	Slide Dryer II, forced air; 220-240VAC, 50/60Hz, 1A; ambient to 75°C; 11"W x 10"D x 5"H



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

[www.trianglebiomedical.com](http://www.trianglebiomedical.com)

Distributed By:

**TBS** TRIANGLE BIOMEDICAL SCIENCES, INC.

3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX: 919.384.9595 • EMAIL: [tbs@trianglebiomedical.com](mailto:tbs@trianglebiomedical.com)

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. TBS is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.

Litho in U.S.A. © March 2007

Bulletin # 12.1.5