

SHUR/Dry™

Make Slide Drying a Breeze



TBS®


SHUR/Dry™ Features:

- Ambient to 80° ± 1°
- 51CFM of air flow
- Digital programming
- Microprocessor controlled
- Fluorescent lighting
- Adjustable shelf
- Four 3,200rpm fans



Four powerful 3,200rpm fans each blow 51cfm of heated air across a 700w, 18” heating strip for rapid drying of slides.

Recessed fluorescent lighting



Magnetic latch


Adjustable shelf

SD-III-120

Removable pan



Adjustable shelf accepts an assortment of slide racks with and/or without coverslipper adapters



SD-II-120

Conventional Dryer Featuring:

- Temperature control to 75° C
- Lexan front for easy viewing
- Thermal switch for overheat protection
- Accommodates two stainer racks

Ordering Information:

Cat #	U/M	Description
SD-II-120	ea	Slide Dryer II, forced air, ambient to 75°C. 100-120VAC, 50/60Hz, 2A; Dimensions WxDxH: (11x10x5)in, (28x26x13)cm
SD-II-220	ea	Slide Dryer II, forced air, ambient to 75°C. 220-240VAC, 50/60Hz, 1A; Dimensions WxDxH: (11x10x5)in, (28x26x13)cm
SD-III-120	ea	SHUR/Dry™ Slide Dryer III, forced air, w/ one adjustable shelf, microprocessor controlled, ambient to 70°C ± 1°, 9W recessed fluorescent light, four 3,200rpm fans. 100-120VAC, 50/60Hz, 2A; External Dimensions WxDxH: (22x16x20)in, (56x41x51)cm. Chamber Dimensions WxDxH: (19x9x14)in, (48x23x36)cm
SD-III-220	ea	SHUR/Dry™ Slide Dryer III, forced air, w/ one adjustable shelf, microprocessor controlled, ambient to 70°C ± 1°, 9W recessed fluorescent light, four 3,200rpm fans. 220-240VAC, 50/60Hz, 1A; External Dimensions WxDxH: (22x16x20)in, (56x41x51)cm. Chamber Dimensions WxDxH: (19x9x14)in, (48x23x36)cm
SD-AS	ea	Shelf for Slide Dryer III, for additional space utilization



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

www.trianglebiomedical.com

Distributed By:

TBS TRIANGLE BIOMEDICAL SCIENCES, INC.

3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX:919.384.9595 • EMAIL: tbs@trianglebiomedical.com

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. TBS is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.