




SHUR/Mark® Cassette Sleeves

*Prepackaged Cassette Sleeves insert directly
into SHUR/Mark® Hoppers*




Innovations for Science™


How To Load SHUR/Mark[®] Cassette Sleeves



Specialty dispensing hopper (SM-TSDH) accepts prepack'd cassette sleeves.




Insert cassette sleeves into hopper sitting on 4-position hopper rack (SM-HR).



Place loaded hopper on SHUR/Mark[®] and pull tape down to seat cassettes.

Ordering Information

	White	Blue	Green	Gray	Pink	Peach	Yellow	Tan	Aqua	Lilac
Cat #	U/M	Description								
SM-PCTS-(color)	cs	SHUR/Mark [®] Process & Embedding cassettes. Prepackaged SHUR/Track [™] taped cassette sleeves for direct insertion into SM-TSDH dispensing hopper receptacle. Colors Available: White (W), Blue (B), Green (G), Gray (GY), Pink (P), Peach (PE), Yellow (Y), Tan (T), Aqua (A) & Lilac (L); Packaged (20 sleeves of 50 cassette bottoms & 1 box of 1000 tops)/cs								
SM-BCTS-(color)	cs	SHUR/Mark [®] Biopsy cassettes. Prepackaged SHUR/Track [™] biopsy taped cassette sleeves for direct insertion into SM-TSDH dispensing hopper receptacle. Colors Available: White (W), Blue (B), Green (G), Gray (GY), Pink (P), Peach (PE), Yellow (Y), Tan (T), Aqua (A) & Lilac (L); Packaged (20 sleeves of 50 cassette bottoms & 1 box of 1000 tops)/cs								
SM-TSDH	ea	Dispensing hopper for SHUR/Track [™] Cassette Tape Sleeves, reusable metal, accepts prepackaged taped cassette sleeves SM-PCTS-(color) and SM-BCTS-(color) biopsy cassettes, 50 cassette capacity								
SM-HR	ea	Hopper Rack, 4 position for holding additional (SM-TSDH) hoppers to facilitate loading of pre-								



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

www.trianglebiomedical.com

Distributed By:

TBS[®] TRIANGLE BIOMEDICAL SCIENCES, INC.
 3014 Croasdaile Drive • Durham, NC 27705-2507 • USA • 919.384.9393 • Fax: 919.384.9595 • E-Mail: tbs@trianglebiomedical.com

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. SHUR/Track[™] is a trademark, and SHUR/Mark[®] & TBS[®] are registered trademarks of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA. Litho in U.S.A. Printed February 2008 Hopper Patent #6.098.839, Patent Applied for Tape Mechanism