


# SHUR/Sharp™

## Disposable Microtome Blade Cutting Solutions



**TBS**®


# The SHUR/Sharp™ Advant“edge” - A cut above the rest



**SHUR/Sharp™** disposable microtome blades are manufactured from premium grade steel for superior blade life and edge retention. State of the art grinding, honing and strapping equipment produce flawless edges. Proprietary coating process minimizes compression and produces excellent ribbons. Select **SHUR/Sharp™** low or high profile blades for consistent sharpness & dependability. Specify **SHUR/Sharp™ PLUS** when ultra thin sections are required.

## Ordering Information:

Cat #	U/M	Description
DMB-LP	pk/cs	<b>SHUR/Sharp™</b> Disposable Microtome Blade, Coated, Low Profile; (0.012" x 0.312" x 3"). 50/pk, 10pk/cs
DMB-LP+	pk/cs	<b>SHUR/Sharp™ PLUS</b> Disposable Microtome Blade, Premium Quality, Ultra Thin, Coated, Low Profile; (0.010" x 0.312" x 3"). 70/pk, 10pk/cs
DMB-HP	pk/cs	<b>SHUR/Sharp™</b> Disposable Microtome Blade, Coated, High Profile; (0.012" x 0.553" x 3"). 50/pk, 10pk/cs
DMB-HP+	pk/cs	<b>SHUR/Sharp™ PLUS</b> Disposable Microtome Blade, Premium Quality Coated, High Profile; (0.012" x 0.553" x 3"). 60/pk, 10pk/cs



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

[www.trianglebiomedical.com](http://www.trianglebiomedical.com)

Distributed By:

**TBS® TRIANGLE BIOMEDICAL SCIENCES, INC.**  
3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX:919.384.9595 EMAIL:tbs@trianglebiomedical.com

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability.

SHUR/Sharp™ is a trademark and TBS® is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.

Litho in U.S.A. Printed February 2006

Bulletin #26,27,4.1