



A Cut Above the Rest



SHUR/Sharp™

A name you can depend on for superior disposable microtome blades, holders and sectioning accessories. The **SHUR/Sharp™** one-step lever activated holder will perform flawlessly for years. Choose a high, low or the new interchangeable high/low profile holder with or without an extended handle to eliminate knife guard interference or select the size that accepts TBS's own proprietary heavy duty blade. We guarantee the same great quality regardless of the blade preference.

SHUR/Sharp™ & SHUR/Sharp™ PLUS

Regular, Heavy Duty, and High & Low Profile Blades

Choose from six disposable microtome blades to meet all your sectioning needs. Standard blades for routine tissue samples and heavy duty blades for cutting difficult tissues and frozens are available for the TBS intermediate height slotted holder. **SHUR/Sharp™** high profile blade fits all high profile holders. **SHUR/Sharp™** high and low profile blades are available to fit any high or low profile blade holder you prefer. **SHUR/Sharp™ PLUS** is the perfect selection when ultra thin sections are required. All blades feature a Teflon coating to reduce compression.



Cryostat Adapters

Allows the Sakura Tissue-Tek® cryostat to take advantage of the benefits of a **SHUR/Sharp™** blade holder. The adapters simply screw onto the back of the blade holder for mounting directly to the microtome. Guaranteed to improve the quality of your sections.

Anti-Roll Device


Snap it on or off a TBS holder. This unique magnetic attachment system positions a clear plastic plate above the blade in order to facilitate specimen handling.



“Our Edge Is The Leading Edge”

Ordering Information

Catalog #	U/M	Description
DMBH	ea	SHUR/Sharp™ Disposable Microtome Bladeholder, 100% Stainless Steel construction, Unconditional 2 year warranty, requires DMH or DMH-HD SHUR/Sharp™ blades.
DMBH-E	ea	Same as DMBH but with extended handle to eliminate knife guard interference.
DMBH-LP	ea	Same as DMBH but designed for use with any low profile microtome blade.
DMBH-LPE	ea	Same as DMBH-LP but with extended handle to eliminate knife guard interference.
DMBH-HP	ea	Same as DMBH but designed for use with DMH-AO or any high profile microtome blade.
DMBH-HPE	ea	Same as DMBH-HP except with extended handle to eliminate knife guard interference.
DMBH-HPLP	ea	Patented SHUR/Sharp™ blade holder designed for interchangeable use of low and high profile blades.
DMBH-CA	ea	Cryostat Adapter required for using DMH (all) in Tissue-Tek® cryostats.
ARD	ea	Anti-Roll Device for use with all TBS® bladeholders.
ARD-P	ea	Replacement plate for ARD.
CA-KS	ea	Knife Slide adapter required for using DMH-HD on the bladeholder provided with Leica Cryostats, models 1800 and 2800.
DMB	pk/cs	SHUR/Sharp™ Disposable Microtome Blade, Coated, Regular (0.012" x 0.498" x 3"). Requires DMH or DMH-E Disposable Microtome Blade Holder. 60/pk, 10pk/cs.
DMB-HD	pk/cs	SHUR/Sharp™ Disposable Microtome Blade, Coated, Heavy Duty, (0.020" x 0.498" x 3"). Requires DMH or DMH-E Disposable Microtome Blade Holder. 35/pk, 10pk/cs.
DMB-LP	pk/cs	SHUR/Sharp™ Disposable Microtome Blade, Coated, Low Profile, (0.012" x 0.312" x 3"). 50/pk, 10pk/cs
DMB-LP+	pk/cs	SHUR/Sharp™ PLUS Disposable Microtome Blade, Coated, Low Profile, (0.010" x 0.312" x 3"). 60/pk, 10pk/cs
DMB-HP	pk/cs	SHUR/Sharp™ Disposable Microtome Blade, Coated, High Profile, (0.012" x 0.553" x 3"). 50/pk, 10pk/cs
DMB-HP+	pk/cs	SHUR/Sharp™ PLUS Disposable Microtome Blade, Coated, Low Profile, (0.012" x 0.553" x 3"). 60/pk, 10pk/cs
TTOH	ea	Moveable Object Holder for Tissue-Tek® cryostat, includes clamp and three (3) ball end chucks.
BEC	ea	Additional ball end chucks for TTOH.



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

www.trianglebiomedical.com

Distributed By:

TBS® TRIANGLE BIOMEDICAL SCIENCES, INC.

3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX: 919.384.9595 EMAIL: tbs@trianglebiomedical.com

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. Tissue Tek® is a registered trademark of Sakura Finetek USA. **SHUR/Sharp™** is a trademark and TBS® is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.

Litho in U.S.A. Printed November 2006

Bulletin #26,27,2.0