

SHUR/Slipper™

One Handed Cover Glass Dispensing



TBS®


Let TBS' SHUR/Slipper™ facilitate your coverslipping

- ◆ Conveniently dispenses cover slips
- ◆ One Hand Operation
- ◆ Portable & Convenient



Ordering Information:

Cat #	U/M	Description
SS-4050	ea/cs	SHUR/Slipper™ Cover Glass Dispensing Instrument for manual coverslipping, accepts prefilled dispensers of SHUR/Slips™ cover glass sizes 24x40 & 24x50; Includes power supply. Specifications: in(mm) 6.5(165) L x 2.4(62) W x 3.4(87) H; 115V, 50/60Hz; 4ea/cs
SS-5560	ea/cs	SHUR/Slipper™ Cover Glass Dispensing Instrument for manual coverslipping, accepts prefilled dispensers of SHUR/Slips™ cover glass sizes 24x55 & 24x60; Includes power supply. Specifications: in(mm) 6.5(165) L x 2.4(62) W x 3.4(87) H; 115V, 50/60Hz; 4ea/cs
SS-2440	pk/cs	SHUR/Slips™ , Prefilled Dispensers of 24x40 cover glass, Thickness No.1 (0.12~0.17mm); approximately 210pc/ea, 5ea/pk, 4pk/cs
SS-2450	pk/cs	SHUR/Slips™ , Prefilled Dispensers of 24x50 cover glass, Thickness No.1 (0.12~0.17mm); approximately 210pc/ea, 5ea/pk, 4pk/cs
SS-2455	pk/cs	SHUR/Slips™ , Prefilled Dispensers of 24x55 cover glass, Thickness No.1 (0.12~0.17mm); approximately 210pc/ea, 5ea/pk, 4pk/cs
SS-2460	pk/cs	SHUR/Slips™ , Prefilled Dispensers of 24x60 cover glass, Thickness No.1 (0.12~0.17mm); approximately 210pc/ea, 5ea/pk, 4pk/cs



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

www.trianglebiomedical.com

Distributed By:

TBS® TRIANGLE BIOMEDICAL SCIENCES, INC.
 3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX:919.384.9595 • EMAIL: tbs@trianglebiomedical.com

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. SHUR/Slipper & SHUR/Slips are trademarks and TBS is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.