

SHUR/Track™ Cassettes

Unique Operation & Unparalleled Security



TBS®

How To Use SHUR/Track™ Cassettes



Insert cover until hinges click firmly into tracks




Cover can remain open once locked into place



Press down on cover & snap securely in place

Ordering Information

White	Blue	Green	Gray	Pink	Peach	Yellow	Tan	Aqua	Lilac
Cat #	U/M	Description							
SM-PPC-(color)	cs	SHUR/Track™ Process & Embedding cassettes w/channel track hinge. Specifically designed and quality controlled for use w/both automated labeling instruments and hand marking devices. Detachable lids are removed and packaged separately. Colors Available: White (W), Blue (B), Green (G), Gray (GY), Pink (P), Peach (PE), Yellow (Y), Tan (T), Aqua (A) & Lilac (L); 1000/cs							
SM-PCTS-(color)	cs	SHUR/Track™ Process & Embedding cassettes. Prepackaged for direct insertion into SM-TSDH dispensing hopper. Colors Available: White (W), Blue (B), Green (G), Gray (GY), Pink (P), Peach (PE), Yellow (Y), Tan (T), Aqua (A) & Lilac (L); Packaged (20 sleeves of 50 cassette bottoms & 1 box of 1000 tops)/cs							
BC-(color)	cs	SHUR/Track™ Biopsy cassettes w/channel track hinge. Specifically designed and quality controlled for use w/both automated labeling instruments and hand marking devices. Detachable lids are removed and packaged separately. Colors Available: White (W), Blue (B), Green (G), Gray (GY), Pink (P), Peach (PE), Yellow (Y), Tan (T), Aqua (A) & Lilac (L); 500/cs							
PPC-BLK-(color)	cs	SHUR/Track™ Process & Embedding Cassettes, same as SM-PPC, except packaged in bulk for large volume applications. Colors Available: White (W), Blue (B), Green (G), Gray (GY), Pink (P), Peach (PE), Yellow (Y), Tan (T), Aqua (A) & Lilac (L); 3000/cs							
PEC-II-(color)	cs	SHUR/Track™ Process & Embedding cassettes, bottoms only, for use w/metal reusable lids. Colors Available: White (W), Blue (B), Green (G), Gray (GY), Pink (P), Peach (PE), Yellow (Y), Tan (T), Aqua (A) & Lilac (L); 500/pk, 3 pk/cs							
SM-BCTS-(color)	cs	SHUR/Track™ Biopsy cassettes w/channel track hinge. Prepackaged for direct insertion into SM-TSDH dispensing hopper. Colors Available: White (W), Blue (B), Green (G), Gray (GY), Pink (P), Peach (PE), Yellow (Y), Tan (T), Aqua (A) & Lilac (L); Packaged (20 sleeves of 50 cassette bottoms & 1 box of 1000 tops)/cs							



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

www.trianglebiomedical.com

Distributed By:

TBS® TRIANGLE BIOMEDICAL SCIENCES, INC.
 3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX:919.384.9595 • EMAIL: tbs@trianglebiomedical.com

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. SHUR/Track™ is a trademark and TBS is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA. Cassette Hinge Patent #6,395,234, Hopper Patent#6,098,839, Patent Applied for Tape Mechanism. Litho in U.S.A. Printed April 2007

Bulletin # 23.3.3