



# SHUR/View™ Microscopes

**Ideal for any Viewing**



*Innovations for Science™*

Whether you need a basic binocular microscope for stain checking and skin margin confirmations or a feature rich instrument for intense DERM and Histopathological applications, TBS can match your requirements with a quality economical solution. Manufactured and assembled in an ISO9001 certified facility, our CE Certified microscope features 22X wide-field eye pieces, infinity plan objectives and too many innovations to list. Contact your TBS representative to determine the features that meet your requirements at an affordable cost.



**SV-3000MOHS**



**SV-3000-IC**



**SV-3000-I**



**Lens Cleaner**

Removes dirt, dust and fogging from lenses while creating a protective layer on lenses to prevent further residue build-up.



**SHUR/View™ Lens Paper**

Dust-free lens tissues are ideal with liquid cleaners. Used to clean microscope lens and objectives.




**SHUR/View™ Immersion Oil**

Type A immersion oil with viscosity of 150 cSt and refractive index 1.5239

**Ordering Information:**

Cat. #	Description
<b>SV-3000-I</b>	<b>SHUR/View™</b> Microscope; Binocular type with 10x22mm eye pieces with caps, infinity plan objectives, 2,4,10,40x, swing out N.A. 1.25 abbe condenser
<b>SV-3000-IC</b>	<b>SHUR/View™</b> Microscope; Trinocular type with 10x22mm eyepieces with caps, infinity plan objectives, 2,4,10,40x swing out 1.25 abbe condenser, 5 Megapixel removable camera.
<b>SV-3000MOHS</b>	<b>SHUR/View™</b> Microscope; Binocular type with 10x22mm eye pieces with caps, 2.5,4,10,20x, infinity plan objectives, swing out 1.25 abbe condenser.
<b>SV-3000MOHS-IC</b>	<b>SHUR/View™</b> Microscope; Trinocular type with 10x22mm eyepieces with caps, with 2.5, 4, 10, 20x infinity plan objectives swing out 1.25 abbe condenser, 5 Megapixel removable digital camera.
<b>SV-LVIO</b>	<b>SHUR/View™</b> Immersion oil, low viscosity; Clear dropper bottle; 1/2oz./ea, 12ea/cs
<b>SV-FIO</b>	<b>SHUR/View™</b> non-Fluorescence Immersion oil; Clear dropper bottle, 1/2oz./ea, 12ea/cs
<b>SV-HVIO</b>	<b>SHUR/View™</b> Immersion oil, high viscosity; Clear dropper bottle, 1/2oz./ea, 12ea/cs
<b>SV-LP</b>	<b>SHUR/View™</b> Lens Paper, 4x6" Booklet, 50 sheets/ea, 5ea/pk
<b>SV-LC</b>	<b>SHUR/View™</b> Lens Cleaner, Spray Bottle, 1oz./ea, 5ea/cs
<b>CUT3400</b>	404 Oil, for sliding surfaces, 100ml.



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

[www.trianglebiomedical.com](http://www.trianglebiomedical.com)

Distributed By:

**TBS® TRIANGLE BIOMEDICAL SCIENCES, INC.**

3014 Croasdaile Drive • Durham, NC 27705-2507 • USA • 919.384.9393 • Fax: 919.384.9595 • E-Mail: [tbs@trianglebiomedical.com](mailto:tbs@trianglebiomedical.com)

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. **SHUR/View™** AND **Innovations for Science™** are trademarks and TBS is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.