

# Purchase Any SHUR/Cut™ Microtome . . .



**SC2500**

## **Semi-Automatic Non-Retracting Microtome**

The SC2500 has the same motorized trimming & sectioning thickness settings features as the SC3500; but, doesn't have the added benefit of an ergonomic hand-pad controller



**SC3500**

## **Semi-Automatic Retracting Manual Microtome**

User-friendly hand-pad allows the operator to store separate trimming and sectioning thicknesses. To switch from one to the other, simply press the appropriate embossed switch. 2-speed motorized specimen arm forward and reverse feed functions rapidly advance a tissue block for cutting or returning to the rear for inserting a new block. Retraction is activated / deactivated by the push of a button. One-handed operation for specimen course advance / return. No longer is the operator required to extend both arms to section or trim. Simply push a button to switch between operator selectable pre-set trimming & sectioning thickness and then rotate the hand-wheel at your preferred speed to generate beautiful ribbons effortlessly.




**SC4500**

## **Automatic Retracting Microtome**

Automatic retracting microtome w/ adjustable-speed motorized cutting stroke. Continuous and single-stroke motorized cutting modes or operate manually for specimens requiring a special touch. Repetitive motion hand functions housed in operator positionable handpad for technician comfort and to reduce work-related Neuro-musculoskeletal disorders. Push-button setting of cutting thickness, switching between trimming/sectioning modes, 2-speed specimen advance/return, O/O retraction feature. Unique specimen position-memory feature invaluable for remembering a specimen's extended position for removing a block to chill and then return automatically to the original position with just a push of a button.

Receive **10% OFF SHUR/Mark®** Cassette Taped Sleeves *OR*

2 cases of **PARA/Gard™** for **FREE!**



### **SM-PCTS-(color) and SM-BCTS-(color)**

**SHUR/Mark®** Processing & Embedding or Biopsy cassette taped sleeves. Prepackaged **SHUR/Track™** cassette sleeves for direct insertion into SM-TSDH taped sleeve dispensing hopper. Colors available: White (W), Blue (B), Green (G), Grey (GY), Pink (P), Peach (PE), Yellow (Y), Tan (T), Aqua (A), & Lilac (L); Packaged (20 sleeves of 50 cassette bottoms & container of 1000 tops)/cs



### **PG-16**

**PARA/Gard™**, liquid paraffin repellent for countertops and equipment; 16 oz. spray bottle/ea, 4 ea/cs

For a downloadable brochures, please visit:

<http://www.trianglebiomedical.com/shurcutseriesmicrotomes.cfm>,

<http://www.trianglebiomedical.com/cassettes.cfm> & <http://www.trianglebiomedical.com/paragard.cfm>

*You must include promo code# 120212Microtome at the time you order your microtome to receive this special limited time offer!*

*Promotion Expires February 29, 2012. Cannot be combined with other promotions.*

