


Purchase a **SHUR/Stain™** Automated Slide Stainer . . .




Comparison Features	TBS SS-3030	Leica Autostainer XL
Number of User-Definable Programs	50	15
Number of Staining Stations	30	18
Simultaneous staining of multiple slide racks	✓	✓
Enclosed system with fume extraction	✓	✓
Number of Drying Stations	2	1
User-Friendly PC based Programming	✓	NO
Separate DI and tap water rinse inlet	✓	NO
Perform Routine & Special Stains	✓	✓
Touchscreen computer for convenient operation	✓	NO
Complete Open-Access to All Wells at One-Time	✓	✓
Compact Design (W • D • H)(in)	34 • 25 • 18	42 • 26 • 20
UL Certification	✓	Unknown
Made in USA 	✓	NO

DISCLAIMER: Based on the information obtained from competitive web-sites & copies of literature, we believe the following to be correct; however, TBS makes no guarantee the information is accurate or current. Please check the latest sources for such information before relaying on the content.

For downloadable brochures, please visit:
<http://www.trianglebiomedical.com/SHURStainStainer.cfm>

Receive up to 10% of your Purchase Price in **FREE SHUR/Stain™** Reagents of your choice!



- AFB Ziehl Neelson Kit (Hot Method)
- Alcian Blue PAS Stain Kit
- Colloidal Iron Stain Kit
- Eosin Y
- Fontana Masson Stain Kit
- Gill's #2 & #3
- Gram's Stain Kit (with Iodine)
- Harris Hematoxylin
- Lithium Blue
- Masson's Trichome Stain Kit
- Prussian Blue Iron Stain Kit
- Periodic Acid Schiff's (PAS) Stain Kit
- Scott's Bluing Reagent
- Verhoeff Van Gieson Elastin Stain Kit



TBS has long been recognized for our industry leading **TMD™** - Tissue Marking Dyes, **SHUR/Mount™** Mounting Mediums, **TFM™** -Tissue Freezing Medium and **POLY/Fin™** Paraffin. We have recently expanded our NC manufacturing facility to include a special stain & bulk chemical packaging operation. You now have available a single source for quality reagents and stains.

*Promo Code# "120212-SS3030" must be referenced at the time of order placement to receive this limited time offer!
 Promotion Expires February 29,2012. Cannot be combined with other special deals.*